

BACTERIA TESTS

DIPSLIDES, TUBE TESTS, INCUBATORS



DIPSLIDES

INCUBATORS

RELIABLE BACTERIAL TESTS

LIQUID AND SURFACE TESTING

FLEXIBLE PADDLE

BACTEDT057
Incubator for 2 dipslides



BACTEDT055
Incubator for 8 dipslides



BACTEDT059
Incubator for 10 dipslides



BACTEDT052
Incubator for 20 dipslides



BACTEDT093
Incubator for 45 dipslides



FEATURES

The definition of a dipslide is a means of testing the microbial content of liquids or surfaces. It consists of a plastic carrier bearing a sterile culture medium which can be dipped in the liquid or pressed on to the surface to be sampled. It is then incubated for 48 hours at a specific temperature to allow microbial growth.

The resulting microbial colonies are estimated by reference to a chart. If the treatment programme is effective, the dipslide count will be consistently low (10^4). If a high result is obtained, the test should be repeated and the treatment programme checked. Water-i.d.® dipslides for bacterial testing are provided as 10 dipslides in a box. Water-i.d.® offers a variety of incubators along with the dipslide range. Each Water-i.d.® Dipslide pack contains 10 dipslides along with a comparison chart to obtain the test result.

DIPSLIDES

ITEM N°	PRODUCT DESCRIPTION
BactD001	TTC/TTC Dipslide Tests for Total Viable Count in waters and surfaces . Flexible paddle for effective surface contact.
BactD002	Malt/TTC Dipslide Tests for Yeasts, Moulds, Fungi and Total Viable Count in waters and surfaces . Flexible paddle for effective surface contact.
BactD003	RBS/TTC Dipslide Tests for Yeasts, Moulds, Fungi and Total Viable Count for food and manufacturing industry . Flexible paddle for effective surface contact.
BactD004	MAC/TTC Dipslide Tests for Coliforms and Total Viable count in process water . Flexible paddle for effective surface contact.
BactD005	E.Coli & Coliforms/TTC Dipslide Tests for E.Coli, Pseudomonas Aeruginosa and Total Viable Count for enumeration of E.Coli . Flexible paddle for effective surface contact.
BactD006	PDM/MAC Dipslide Tests for Pseudomonas – Coliforms/Enterobacterceae for potable water, pools and spas . Flexible paddle for effective surface contact.
BactD007	PDM/TTC Dipslide Tests for Pseudomonas species and Total Viable Count for Closed water systems, pools and spas . Flexible paddle for effective surface contact.
BactD008	SRB Tube Tests for Sulphide reducing bacteria in Closed water systems . Flexible paddle for effective surface contact.
BactD009	NRB Tube Tests for Nitrite reducing bacteria in Closed water systems . Flexible paddle for effective surface contact.
BactD011	BactD011 R2A/R2A Dipslide Tests for A. Hydrophila, P. Fluorescens, Strep. Pyogenes and E.Coli for potable water and waste water

INCUBATORS

ITEM CODE	PRODUCT DESCRIPTION	VOLTAGE
BactEDT057	Incubator for 2 Dipslide Tests . Preset to 30°C.	10 220 V
BactEDT055	Incubator for 8 Dipslide Tests . With temperature gauge and adjustable temperature.	12 110 220 V
BactEDT059	Incubator for 10 Dipslide Tests . With temperature gauge and adjustable temperature.	12 110 220 V Ideal for NRB/SRB tests. Only offered with 240 V.
BactEDT052	Incubator for 20 Dipslide Tests . With temperature gauge and adjustable temperature.	12 110 220 V
BactEDT093	Incubator for 45 Dipslide Tests . With temperature gauge and adjustable temperature.	12 110 220 V
BactEDT-12V	12V car-lead for BactEDT057 BactEDT055 BactEDT052 BactEDT093	12 V



DIPSLIDE PACKS

Dipslide Packs are always 10 double sided dipslides with flexible paddle for effective surface contact.
Package includes comparison chart and instructions.

Item No.	Product Description
BactD001	TTC/TTC Dipslide Tests for Total Viable Count in waters and surfaces.
BactD002	Malt/TTC Dipslide Tests for Yeasts, Moulds, Fungi and Total Viable Count in waters and surfaces.
BactD003	RBS/TTC Dipslide Tests for Yeasts, Moulds, Fungi and Total Viable Count for food and manufacturing industry.
BactD004	MAC/TTC Dipslide Tests for Coliforms and Total Viable count in process water.
BactD005	Chromogenic/TTC Dipslide Tests for E.coli, Pseudomonas Aeruginosa and Total Viable Count for enumeration of E.coli.
BactD006	PDM/MAC Dipslide Tests for Pseudomonas species and Coliforms for potable water, pools and spas.
BactD007	PDM/TTC Dipslide Tests for Pseudomonas species and Total Viable Count for Closed water systems, pools and spas.
BactD008	SRB Tube Tests for Sulphate reducing bacteria in Closed water systems.
BactD009	NRB Tube Tests for Nitrite reducing bacteria in Closed water systems.
BactD011	R2A / R2A Dipslide Tests for A. Hydrophila, P. Fluorescens, Strep. Pyogenes and E.Coli for potable water and waste water



INCUBATORS

Item No.	Product Description	Maximum dipslide qty	Voltage
BactEDT057	Incubator for 2 Dipslide Tests. Preset to 30°C.	2	110 / 220 V
BactEDT055	Incubator for 8 Dipslide Tests. Preset to 30°C.	8	12 / 110 / 220 V
BactEDT059	Incubator for 10 Dipslide Tests. With temperature gauge and adjustable temperature. Ideal for NRB/SRB tests.	10	12 / 110 / 220 V
BactEDT052	Incubator for 20 Dipslide Tests. With temperature gauge and adjustable temperature.	25	12 / 110 / 220 V
BactEDT093	Incubator for 45 Dipslide Tests. With temperature gauge and adjustable temperature.	45	12 / 110 / 220 V
BactEDT-12V	12V car-lead for BactEDT057/BactEDT055/ BactEDT052/BactEDT093	--	12 V

We gladly Private Label our products for you! Write to sales@water-id.com to receive your individual quotation